



L-17 Large Magnetic Stirrer 60 ltr.

TECHNICAL DATA

- Topplate 320 x 330 mm
Stainless steel 50-1000 rpm
- Dimensions (WxDxH) 320 x 340 x 105 mm
- Art.nr.17000 (230V) 17010 (115V)
- Stirring capacity 60 Ltr.
- Weight 5.2 kg
- incl. spinbar

• Number of positions	1
• Heating quantity max	60 ltr
• Motor rating input(W)	50 W
• Motor rating output(W)	20 W
• Speed display	analog
• Speed range(rpm)	50 - 1000 rpm
• Stirring bar length max(mm)	120 mm
• Heat output	-
• Heating rate	-
• Heating temperature range	-
• Heat control	-
• Heat control accuracy	-
• Speed control	stepless
• Adjustable safety circuit min	-
• Adjustable safety circuit max	-

• Connection for ext. temperature sensor	-
• Control accuracy with sensor	-
• Temperature constancy in medium	-
• Set-up plate material	stainless steel
• Set-up plate dimensions	320 x 330 mm
• Dimensions (W x D x H)	320 x 340 x 105 mm
• Weight(kg)	5.2 kg
• Permissible ambient temperature	5 - 40 °C
• Permissible relative humidity	80%
• Protection class according to DIN EN 60529	IP 42
• RS 232 interface	-
• Analog output	-
• Voltage	220 - 230 V (115 / 100 V)
• Frequency	50/60 Hz
• Power input	50 W

• See also the other models for more options on www.labinco-bv.com



LABINCO B.V.

P.O. BOX 3259 • 4800 DG Breda, The Netherlands
Zinkstraat 64 • 4823 AC Breda, The Netherlands
T 0031(0)76 - 5417582 E info@labinco-bv.com
www.labinco-bv.com