

230v

Multi-position
Design

Vibro shakers

Single-position

Quality

L-26 Small Rotary Mixer

Manufacturing in the Netherlands

Peristaltic pumps

Aluminium topplate

Laboratory equipment

Ceramic coating

Micro-plate shakers

TECHNICAL DATA

Heating-blocks

Rotary mixers

Hotplate Magnetic Stirrers

- Small Rotary Mixer for serial rotation of test-tubes or drums
- For bars or drums see www.labinco-bv.com
- Quick coupling for exchange of rotating drum or bar
- Dimensions (WxDxH) 540 x 160 x 185 mm
- Freq., power input
- Weight 7.5 kg
- Speed 1-50 rpm, analog
- Provided with stainless steel drip tray
- Art.nr.26000 (230V) 26010 (115V)

• Number of positions	1 bar or drum holder
• Mixing quantity max	8 up to 27 mm or max 2 kg
• Motor rating input(W)	40 W
• Motor rating output(W)	25 W
• Speed display	analog
• Speed range(rpm)	1 - 50 rpm
• Stirring bar length max(mm)	-
• Heat output	-
• Heating rate	-
• Heating temperature range	-
• Heat control	-
• Heat control accuracy	-
• Speed control	stepless
• Adjustable safety circuit min	-

• Adjustable safety circuit max	-
• Connection for ext. temperature sensor	-
• Control accuracy with sensor	-
• Temperature constancy in medium	-
• Set-up plate material	steel / nylon
• Set-up spring clip dimensions	8 up to 27 mm
• Dimensions (W x D x H)	540 x 160 x 185 mm
• Weight(kg)	7.5 kg
• Permissible ambient temperature	5 - 40 °C
• Permissible relative humidity	80%
• Protection class according to DIN EN 60529	IP 42
• RS 232 interface	-
• Analog output	-
• Voltage	220 - 230 V (115 / 100 V)

- See also the other models for more options on www.labinco-bv.com



LABINCO B.V.

P.O. BOX 3259 • 4800 DG Breda, The Netherlands
 Zinkstraat 64 • 4823 AC Breda, The Netherlands
 T 0031(0)76 - 5417582 E info@labinco-bv.com
www.labinco-bv.com