



L-32 Basic Hotplate Magnetic Stirrer

20 ltr. 325°C

TECHNICAL DATA

- Hotplate 195 mm diam. Aluminium 50-1250 rpm
- Hotplate temperature 50-325°C
- Dimensions (WxDxH) 200 x 250 x 135 mm
- Art.nr.32000 (230V) 32010 (115V)
- Stirring capacity 20 Ltr.
- Weight 4.2 kg
- incl. spinbar

• Number of positions	1
• Heating quantity max	20 ltr
• Motor rating input(W)	35 W
• Motor rating output(W)	15 W
• Speed display	analog
• Speed range(rpm)	50 - 1250 rpm
• Stirring bar length max(mm)	80 mm
• Heat output	600 W
• Heating rate	7 K/min
• Heating temperature range	room temp. -325 °C
• Heat control	analog
• Heat control accuracy	1 ±K
• Speed control	analog
• Adjustable safety circuit min	-
• Adjustable safety circuit max	-

• Connection for ext. temperature sensor	-
• Control accuracy with sensor	-
• Temperature constancy in medium	1 ±K
• Set-up plate material	aluminium
• Set-up plate dimensions	Ø 195 mm
• Dimensions (W x D x H)	200 x 250 x 135 mm
• Weight(kg)	4.2 kg
• Permissible ambient temperature	5 - 40 °C
• Permissible relative humidity	80%
• Protection class according to DIN EN 60529	IP 42
• RS 232 interface	-
• Analog output	-
• Voltage	220 - 230 V (115 / 100 V)
• Frequency	50/60
• Power input	635 W

• See also the other models for more options on www.labinco-bv.com



LABINCO B.V.

P.O. BOX 3259 • 4800 DG Breda, The Netherlands
Zinkstraat 64 • 4823 AC Breda, The Netherlands
T 0031(0)76 - 5417582 E info@labinco-bv.com
www.labinco-bv.com