



L-366 6-Pos. Flask Heater 450°C

TECHNICAL DATA

- 6 x hotplate 85 mm diam. steel coated
- Hotplate temperature 50-450°C
- Dimensions (WxDxH) 900 x 240 x 130 mm
- Aluminium flask holders
(not included see www.labinco-bv.com)

- Weight 11.2 kg
- Art.nr.36600 (230V)

• Number of positions	6
• Heating quantity max	100-1000 ml (see flask holders)
• Motor rating input(W)	-
• Motor rating output(W)	-
• Speed display	-
• Speed range(rpm)	-
• Stirring bar length max.	-
• Heat output	450 W (1 position)
• Heating rate	7 K/min
• Heating temperature range	room temp. - 450 °C
• Heat control	analog
• Heat control accuracy	1 ±K
• Speed control	-
• Adjustable safety circuit min	-
• Adjustable safety circuit max	-

• Connection for ext. temperature sensor	-
• Control accuracy with sensor	-
• Temperature constancy in medium	1 ±K
• Set-up plate material	steel coated
• Set-up plate dimensions	Ø 85 mm
• Dimensions (W x D x H)	900 x 240 x 130 mm
• Weight(kg)	11.2 kg
• Permissible ambient temperature	5 - 60 °C
• Permissible relative humidity	80%
• Protection class according to DIN EN 60529	IP 42
• RS 232 interface	-
• Analog output	-
• Voltage	220 - 230 V
• Frequency	50/60 Hz
• Power input	2700 W

• See also the other models for more options on www.labinco-bv.com



LABINCO B.V.

P.O. BOX 3259 • 4800 DG Breda, The Netherlands
Zinkstraat 64 • 4823 AC Breda, The Netherlands
T 0031(0)76 - 5417582 E info@labinco-bv.com
www.labinco-bv.com