



**L-46**

## Power Mixer

### TECHNICAL DATA

- Strong vibrating action, for test tubes 7 mm diam. up to 35 mm
- Press-to-Mix action and continuous action
- Dimensions (WxDxH) 135 x 120 x 135 mm
- 4 mm Amplitude
- Weight 3.2 kg
- Art.nr.46000 (230V) 46010 (115V)
- See additional accessorier for a complete mixer line from one basis unit on [www.labinco-bv.com](http://www.labinco-bv.com)

|                                 |                                  |
|---------------------------------|----------------------------------|
| • Number of vibrating positions | 1                                |
| • Vibrating quantity max        | test tubes 7 mm up to 35 mm Dia. |
| • Motor rating input(W)         | 40 W                             |
| • Motor rating output(W)        | 25 W                             |
| • Speed display                 | analog                           |
| • Speed range(rpm)              | 50 - 1250 rpm                    |
| • Stirring bar length max.      | -                                |
| • Heat output                   | -                                |
| • Heating rate                  | -                                |
| • Heating temperature range     | -                                |
| • Heat control                  | -                                |
| • Heat control accuracy         | -                                |
| • Speed control                 | stepless                         |
| • Adjustable safety circuit min | -                                |
| • Adjustable safety circuit max | -                                |

|                                              |                           |
|----------------------------------------------|---------------------------|
| • Connection for ext. temperature sensor     | -                         |
| • Control accuracy with sensor               | -                         |
| • Temperature constancy in medium            | -                         |
| • Set-up plate material                      | rubber cup                |
| • Set-up plate dimensions                    | 7 mm up to 35 mm Diam.    |
| • Dimensions (W x D x H)                     | 135 x 120 x 135 mm        |
| • Weight(kg)                                 | 3.2 kg                    |
| • Permissible ambient temperature            | 5 - 40 °C                 |
| • Permissible relative humidity              | 80%                       |
| • Protection class according to DIN EN 60529 | IP 42                     |
| • RS 232 interface                           | -                         |
| • Analog output                              | -                         |
| • Voltage                                    | 220 - 230 V (115 / 100 V) |
| • Frequency                                  | 50/60 Hz                  |
| • Power input                                | 50 W                      |

See also the other models for more options on [www.labinco-bv.com](http://www.labinco-bv.com)



#### LABINCO B.V.

P.O. BOX 3259 • 4800 DG Breda, The Netherlands  
 Zinkstraat 64 • 4823 AC Breda, The Netherlands  
 T 0031(0)76 - 5417582 E [info@labinco-bv.com](mailto:info@labinco-bv.com)  
[www.labinco-bv.com](http://www.labinco-bv.com)