

Magnetic Stirrers

Model	Model L-71	Model L-73
Cat. number	71000	73000
Description	Economy Magnetic Stirrer	Large Capacity Magnetic Stirrer
Stirring surface	155 mm diam. Stainless Steel	
Stirring capacity	5 ltr.	20 ltr.
Motor rating	22/2,2 W	13/5 W
Speed range	100-1250 rpm	
Weight	2,2 kg	
Dimensions (wx dxh)	172 x 286 x 80 mm	
Permissible ambient	5 - 50° C	
Permissible humidity	90%	
Protection class	IP 42	
Voltage/Frequency	220V/50Hz	



Model L23 Magnetic-Stirrer

- Electronic speed control
- Speed; 50 - 1100 rpm
- Housing: 2 piece casted aluminum
- Dripwater-spillage proved design and construction
- Electrostatically applied powder coated
- M-10 threaded support rod connection
- Supplied with teflon magnetic stirring bar line cord and manual
- 10 liter capacity stainless steel topplate: 180 mm ø

Model	Cat. number	dimensions (l x d x h)	Weight
LD-23	23000 (230V)	200 x 250 x 90 mm	2,9 kg
	23010 (115V)	200 x 250 x 90 mm	2,9 kg



LABINCO B.V.

P.O. BOX 3259 - 4800 DG Breda, The Netherlands

Zinkstraat 64 - 4823 AC Breda, The Netherlands

Tel.: 0031(0)76 - 5417582, Fax: 0031(0)76 - 5420173

E-mail: info@labinco-bv.com - Internet: www.labinco-bv.com