

Magnetic Stirrers



Model L22 LUM-Stirrer

- Opal polycarbonate topplate, 160 mm, diam.
- Two gas-discharge tubes for liquid illumination
- Electronic speed control
- M-10 threaded support ROD connection
- 230 V 50 Hz only

Available	Cat. number
Single-position	22000
4 position - unit	22400



Large Stirrer 35 liter

- Optical speed-detection
- Reversible stirring
- Topplate st.st.
- Micro-processor controlled

Model	Cat. number	mode	dimensions (lxdxh)	Weight
LD-14	14000	digital	320 x 340 x 105 mm	5,2 kg
L-17	17000	not digital	320 x 340 x 105 mm	5,2 kg