

Rock'n Rollers

Available in two models:

- **Single Speed Model L-201**
- **Single Regulated Model L-202**

Both models provide:

- Gentle rocking while rolling, to keep samples in motion
- Rocking movement 35 degr.
- 14 stainless steel rotating bars 208 x 10 mm diam.
- Heart to heart distance of tubes 19,5 mm



Single-speed Rock'n roller Model L-201

- Roller speed: 10 rpm
- Rocking-frequency: 5 Cy/min
- Dimensions: 320 x 310 x 190 mm
- Weight: 3,6 kg

Ordering:

230V 50/60Hz. Cat.no. 20100
115 V 50/60Hz. Cat.no. 20110



Variable-speed Rock'n roller Model L-202

- Rolling Speed: from 5 to 25 rpm
- Extension Cable: 2 m for in-cabinet use
- Rocking-frequency: 2,5 to 12,5 Cy/min
- Temperature Range:
in cooler or incubator 2 - 50°C
- Basic Unit:
Dimensions: 320 x 310 x 190 mm
Weight 3,6 kg
- Controller:
Dimensions: 150 x 190 x 75 mm
Weight 1,2 kg

Ordering:

230V 50/60Hz. Cat.no. 20200
115 V 50/60Hz. Cat.no. 20210



LABINCO B.V.

P.O. BOX 3259 - 4800 DG Breda, The Netherlands

Zinkstraat 64 - 4823 AC Breda, The Netherlands

Tel.: 0031(0)76 - 5417582, Fax: 0031(0)76 - 5420173

E-mail: info@labinco-bv.com - Internet: www.labinco-bv.com