

Multi-position-magnetic-stirring-systems



The "Labinco" Models L-714/744/716/746/ is a 4 or 6 position instrument that provides stirring in beaker glasses, erlenmeyers or other suitable containers.

- Stainless steel toplate
- Stirring speed is stepless controlled between 0-1200 rpm (clockwise/anti-clockwise)

Specifications; 230 Volt	4 pos. =100 Watt	6 pos. =150 Watt	50/60Hz	(other lines on request)
BENCH-MODELS L-714 / 716				
Dimensions ; (L x W x H)	L-714= 545 x 245 x 155 mm (weight 8 kg)			cat.no: 71400
	L-716= 845 x 245 x 155 mm (weight 12 Kg)			cat.no: 71600
SQUARE-MODELS L-744 / 746				
	L-744= 355 x 400 x 255 mm (weight 8 Kg)			cat.no: 74400
	L-746= 505 x 400 x 255 mm (weight 12Kg)			cat.no: 74600



LABINCO B.V.

P.O. BOX 3259 - 4800 DG Breda, The Netherlands

Zinkstraat 64 - 4823 AC Breda, The Netherlands

Tel.: 0031(0)76 - 5417582, Fax: 0031(0)76 - 5420173

E-mail: info@labinco-bv.com - Internet: www.labinco-bv.com