

Multi-position-hotplate-magnetic-stirring-systems



The "LABINCO" Models LD-814/844/816/846 is a 4 or 6 position instrument that provides separate or simultaneous heating and stirring in beaker glasses, erlenmeyers or other suitable containers.

- Stainless steel hotplates $\varnothing$ 120 mm
- With digital temperature controller up to 420 C°
- Stirring speed is stepless controlled between 50-1200 rpm



Specification; 230 Volt	4 pos. = 1300 watt	6 pos. =1950 watt	50/60 Hz
	BENCH-MODELS Ld-814/816		
Dimensions; (L x W x H)	Ld-814 = 545 x 245 x 155 mm (weight 10 kg)	cat.no:81400	
	Ld-816 = 845 x 245 x 155 mm (weight 14 kg)	cat.no:81600	
	SQUARE-MODELS LD-844/846		
	Ld-844 = 355 x 400 x 255 (weight 10 kg)	cat.no:84400	
	Ld-846 = 505 x 400 x 255 (weight 14 Kg)	cat.no:84600	