




LD-76 Digital Rotary Mixer

TECHNICAL DATA

- Digital Rotary Mixer for small serial rotation of test-tubes 4-27 mm
- For bars or drums see www.labinco-bv.com
- Quick coupling for exchange of rotating drum or bar
- Dimensions (WxDxH) 490 x 160 x 205 mm
- Weight 8 kg
- Speed 1-99 rpm, digital
- With timer
- Memory function and Re-start function
- Art.nr.76000 (230V) 76010 (115V)

• Number of positions	1 bar or drum holder
• Mixing quantity max	1 kg
• Motor rating input(W)	50 W
• Motor rating output(W)	35 W
• Speed display	digital
• Speed range(rpm)	1 - 99 rpm
• Timer display	digital 1s. up to 99h 59m 59s
• Heat output	-
• Heating rate	-
• Heating temperature range	-
• Heat control	-
• Heat control accuracy	-
• Speed control	digital
• Adjustable safety circuit min	-
• Adjustable safety circuit max	-

• Connection for ext. temperature sensor	-
• Control accuracy with sensor	-
• Temperature constancy in medium	-
• Set-up plate material	steel / nylon clamps
• Set-up spring clip dimensions	see www.labinco-bv.com
• Dimensions (W x D x H)	540 x 160 x 235 mm
• Weight(kg)	8 kg
• Permissible ambient temperature	5 - 40 °C
• Permissible relative humidity	80%
• Protection class according to DIN EN 60529	IP 42
• RS 232 interface	-
• Analog output	-
• Voltage	220 - 230 V (115 / 100 V)
Frequency	50/60 Hz
Power input	60 W

• See also the other models for more options on www.labinco-bv.com



LABINCO B.V.

P.O. BOX 3259 • 4800 DG Breda, The Netherlands
Zinkstraat 64 • 4823 AC Breda, The Netherlands
T 0031(0)76 - 5417582 E info@labinco-bv.com
www.labinco-bv.com